

3M Infection Prevention Solutions
3M™ Comply™ Lead Free Indicator Tape
for Steam Sterilization

Not Made with Natural Rubber Latex



Lead Free

Solvent Free Hassle Free

Lead
the Way

Why it is so important to “get the lead out.”

Standardizing on lead-free products, such as 3M™ Comply™ Steam Indicator Tapes, makes good sense – by helping to simplify operations, avoid the risk of penalties for non-compliance, and enhancing your reputation for good environmental stewardship.

3M adhesives are engineered to fit our customers' needs and are used in such diverse products as airplanes and cell phones. The industry-standard product for exposure monitoring now has even more to offer.

- Does not require hazardous waste disposal
- Redesigned packaging uses up to 30% less plastic
- Convenient, individually wrapped, easy-to-open rolls
- Sustainable, solvent-free manufacturing process



Ordering Information

Indicator Tapes

These Comply Indicator Tapes (ANSI/AAMI/ISO Class 1 Process Indicators) show *at-a-glance* whether packs have been exposed to a steam sterilization process and function as a closure to securely seal the sterilization packs. They accept hand-writing and labeling for recordkeeping information.

Catalog Number	Product Name	Size	Rolls/Intermediate	Rolls/Case
1322-12MM	3M™ Comply™ Lead Free Steam Indicator Tape (designed for use on all wraps)	0.47 in x 60 yd (12mm x 55m)	–	42
1322-18MM		0.70 in x 60 yd (18mm x 55m)	–	28
1322-24MM		0.94 in x 60 yd (24mm x 55m)	–	20
1322-48MM		1.89 in x 60 yd (48mm x 55m)	–	10
1355-18MM	3M™ Comply™ Lead Free Steam Indicator Tape (designed for use on disposable wraps)	0.70 in x 60 yd (18mm x 55m)	–	28
1355-24MM		0.94 in x 60 yd (24mm x 55m)	–	20

Indicator Tape Dispenser

Catalog Number	Product Name	Description	Items/Box	Boxes/Case
C22	Heavy Duty Tape Dispenser	Holds 2 rolls	1	2
M52	Heavy Duty Tape Dispenser	With Tabber	1	1


Best Practices

Always place an AAMI Class 1 Process Indicator* on the outside of every pack, wrapped tray, container or peel pouch, unless the internal chemical indicator is readily visible.

Unexposed




Exposed



Unexposed



Exposed



The chemical indicator ink, printed in a diagonal striped pattern, changes color or turns darker when exposed to steam.

*As defined by ANSI/AAMI/ISO11140-1:2005 (R)2010

Sustainable Solutions for a Healthier World

Environmental Stewardship

In partnership with communities around the globe, we continually strive to protect our planet's health by developing sustainable solutions, business practices and manufacturing processes. Protecting the environment is central to the mission of 3M.

Economic Success

By creating viable business opportunities, investing in local economies and supporting growth and development in emerging markets, we're always looking for innovative ways to enrich the vitality of local communities around the world.

Social Responsibility

We are committed to fostering vibrant communities by actively contributing our expertise, resources and financial support on both corporate and individual levels.

Polyscientific Enterprise Sdn Bhd

百利企业有限公司 152064 D

Melaka (HQ)
272, Taman Asean,
Jalan Malim,
75250 Melaka

T +606 3350690 F +606 3351631

Shah Alam
14, Jalan Serendah 26/39,
Kawasan Perindustrian HICOM,
40400 Shah Alam, Selangor

T +603 51036920 F +603 51036980

Penang
16, Jalan Sungai Tiram 7,
11900 Bayan Lepas,
Penang

T +604 6371500 F +604 6371600

Johor
No. 365 Jalan Ekoperniagaan 11,
Taman Ekoperniagaan 2,
Senai Airport City, 81400 Senai,
Johor

T +607 5955243 F +607 5955745